



BAY-GEN

Institute of Plant Genomics, Human Biotechnology and Bioenergy



Capillary Sequencing Platform



CAPILLARY SEQUENCING

- You provide DNA sample with a purity of 1.8-2.0 (OD260/280)
PCR fragment: $\geq 20 \mu\text{L}$, 50 ng/ μL
plasmid: $\geq 20 \mu\text{L}$, 100 ng/ μL
- We perform cycle sequencing with ABI BigDye[®] Terminator v3.1 using custom or universal primers
removal of unincorporated nucleotides
capillary electrophoresis using ABI3500 Genome Analyzer
- We guarantee >500 base high-quality read length (standard plasmid)
within 3 working days



NET PRICES

1 sample	7.5 Euro
100 samples (pre-paid package)	6.9 Euro/sample
500 samples (pre-paid package)	6.7 Euro/sample
1000 samples (pre-paid package)	6.3 Euro/sample

Send your order and samples:
Capillary Sequencing Platform, BAYGEN
6726 Szeged, Derkovits fasor 2., Hungary
Tel.: +36 30 827 0457
Fax: +36 62 546 974
www.baygen.hu
seq@baygen.hu

